


Biological Indicator



The Biological Indicator Gasket

Rubber Fab is proud to introduce a new product that enables the use of a **self contained biological indicator** in sanitary process lines. This unique Rubber Fab Smart Gasket® permits either top loading or inline positioning of an **EZTest® self contained biological indicator and/or a temperature probe**.

The great advantage of a self contained biological indicator over the traditional spore strip in glassine is the time saved when confirming spore kill. ***Results are easily obtained in only 24 hours (versus the 5 to 7 days required with a spore strip).***

You can feel secure knowing that the biological indicator will not be lost down stream in the process. The fail-safe Rubber Fab Smart Gasket® holds the indicator in place for easy retrieval. Available for 3/4", 1.0", 1.5" and 2.0" Sanitary, Tri-Clamp® process lines.

Meets High Pharmaceutical Standards

- U.S. Pharmacopeia Class VI Certification
- Cytotoxicity Criteria
- CFR Title 21 Section 177.2600
- Traceability: Lot and Batch
- Certification: Lot and Batch
- ASME-BPE Standards
- USDA Standards
- 3-A Sanitary Standards
- Current Good Manufacturing Practices (CGMP)
- Metal Detectable Elastomers available
- Animal Derived Ingredient Free

EZTest®

- U.S. Pharmacopeia
- ISO 11138
- European Norm 866